

# Intra Uterine Contraception



01

All copper IUDs may be used as Emergency Contraceptives, inserted up to 5 days after unprotected intercourse.

All copper IUDs fitted to a woman over the age of 40 years may be left in situ until after the menopause. Recommended replacement intervals shown refer to the use of IUDs in younger women.

## M Shape

### Load 375 01

Copy of the familiar 'M-Shaped' device with 375mm<sup>2</sup> copper, for uterine length over 7cm.

<b>Replace every 5 years</b>	<b>IUD0073001</b>	
1 to 12		each £7.95
13 to 49		each £7.56
50+		each £7.40



02

## Frameless

### GyneFix Original 02


Unique frameless device with 330mm<sup>2</sup> copper as tubes suspended on a surgical suture, the knot of which is embedded into the myometrium – special training required for fitting. Suitable for all uterine lengths.

<b>Replace every 5 years</b>	<b>IUD0041200</b>	<b>£26.64</b>
------------------------------	-------------------	---------------

### GyneFix 200

Features just 4 copper beads providing 200mm<sup>2</sup> copper.

<b>Replace every 5 years</b>	<b>IUD0041300</b>	<b>£26.64</b>
------------------------------	-------------------	---------------



03

## Plain Copper T

### UT 380 Standard 03

Copy of the familiar 'copper T' design, loaded by pulling on the strings to draw the arms back into the introducer, 380mm<sup>2</sup> copper on vertical stem only. For uterine lengths 6.5 to 9cm.

<b>Replace every 5 years</b>	<b>IUD0072001</b>	
1 to 12		each £8.95
13 to 49		each £8.56
50+		each £8.40

### UT 380 Short

As above, but now even smaller. Provides 380mm<sup>2</sup> copper suitable for minimum uterine length of 5cm.

<b>Replace every 5 years</b>	<b>IUD0072002</b>	
1 to 12		each £8.95
13 to 49		each £8.56
50+		each £8.40

# Intra Uterine Contraception

Copper T 380 devices with additional bands of copper on the horizontal arms give the highest efficacy and longest life in the patient. That's why they are regarded as 'The Gold Standard' by the leading authorities<sup>1,2,3</sup> and meet the NICE guidelines for LARC.<sup>4</sup>

## Banded Copper T

### TT 380 Slimline 04

'Banded copper T' design – the 'gold standard', featuring copper on horizontal arms as well as on the vertical stem, 380mm<sup>2</sup> copper, with easy-loading system, no capsule. For uterine length 6.5 to 9cm.

**Replace every 10 years**

**IUD0071001**

**1 to 12**

**each £9.95**

**13 to 49**

**each £9.56**

**50+**

**each £9.40**



04

### Mini TT 380 Slimline

Mini TT 380 Slimline – the only banded copper T device designed for use in the nulliparous. Features copper sleeves on horizontal arms as well as on the vertical stem, providing 380mm<sup>2</sup> copper on a short stem and with narrower arms. Suitable for minimum uterine length of 5cm.

**Replace every 5 years**

**IUD0074001**

**1 to 12**

**each £9.95**

**13 to 49**

**each £9.56**

**50+**

**each £9.40**



04

A demonstration video of how to load the TT380 Slimline is available on our website.

## Intra Uterine System

### Mirena 05

Controlled release of levonorgestrel directly into the uterus for 5 years. Also licensed for menorrhagia and endometrial protection. (N.B. Mirena is not effective as an emergency contraceptive). Replace every 5 years if prescribed for contraception or menorrhagia and every 4 years if prescribed for endometrial protection. Recommended by NICE as one of the LARC options.

**IUS0011100**

**£85.66**



05

#### References:

1. Professor J Guillebaud, Reviews in Gynaecological Practice, 1, (1), 2001, 24-31
2. Jo Dennis & Naomi Hampton, The Journal of Family Planning & Reproductive Healthcare, 28, (2), 2002, 61-68
3. Sunanda Gupta & Rosemary Kirkman, The Obstetrician & Gynaecologist, 1, (4), 2002, 37-43
4. NICE Clinical Guideline 30, October 2005